

STERYLAB PRODUCT RANGE

*Biopsy & localisation needles
Gynaecological products
Decontaminating mats for the O.R.
Aghi per biopsia e localizzazione
Prodotti per ginecologia
Tappeti decontaminanti*

For further information or a quotation please contact : Colin Walters



SeeDOS[™]

Radiation Oncology, Urology and General Medical Equipment and Services

26, The Maltings, Leighton Buzzard, Bedfordshire LU7 4BS, United Kingdom

Tel: +44 1525 850670 Fax: +44 1525 850685 Email: cwalters@seedos.com

Website: www.seedos.com Mobile: +44 07701059163

**DISPOSITIVI PER BIOPSIA DELLA PELLE
SKIN BIOPSY PUNCHES & CURETTES**

DERM PUNCH[®]



Lama in acciaio inossidabile
con affilatura a rasoio

■ *Stainless steel sharp
cutting edge*




■ Indicazione diametro lama

■ *Size embossed on every punch*



DISPOSITIVI PER BIOPSIA DELLA PELLE

SKIN BIOPSY PUNCHES & CURETTES


 La biopsia della pelle, ormai procedura essenziale per la diagnosi istologica di:

- Dermatosi dubbie
- Lesioni sospette di tumore
- Nei anomali

nonché nella diagnosi di infezioni batteriche o micotiche, stati infiammatori patologici, di alcune specifiche malattie come la psoriasi, per verificare l'assenza di cellule tumorali in zona di rimozione del tumore, per monitorare l'efficacia di un trattamento, per rimuovere porri, nei o anomale escrescenze, in alcune lesioni della mucosa orale.

► Dermo-Punch®

- Lama in acciaio inossidabile estremamente affilata
- Impugnatura anti-scivolo per una sicura presa e controllo
- Diametro della lama stampato in rilievo sull'asta per una sicura identificazione
- Vasta gamma dimensionale (tabella)
- Disponibile sia come Punch che Curetta
- Sterile, monouso: riduzione del rischio di infezioni nosocomiali.

 *Skin biopsy, a very common procedure to date, is essential for histologic diagnosis of:*

- *Obscure dermatoses*
- *Lesions suggesting malignancy*
- *Suspected moles*

and to diagnose bacterial or fungal infection, inflammatory skin disorders (such as psoriasis), benign skin grow, to verify the presence of normal skin at the edge of an area from which a tumor was removed, to monitor the effectiveness of a treatment, to remove warts, moles, skin cancers, or other growths and some oral mucosal lesions as well.

► Dermo-Punch®

- *Stainless steel razor sharp cutting edge*
- *Ribbed handle for sure grip and control*
- *Size embossed on each punch for easy identification*
- *Wide range of sizes*
- *Available either as Punch and Curette*
- *Sterile disposable: reduction of infection hazard.*

MODELLI DISPONIBILI

AVAILABLE SIZES

DERMO PUNCH®

TIPO Type	CODICE Cat. Nr.	DIAMETRO LAMA BLADE SIZE (mm)	CONFEZ. PCKG pz-units
PUNCH	DMP02	2	25
	DMP03	3	
	DMP35	3,5	
	DMP04	4	
	DMP05	5	
	DMP06	6	
	DMP08	8	
CURETTE	DMP02C	2	
	DMP03C	3	
	DMP04C	4	
	DMP05C	5	
	DMP07C	7	

Emiss.: D.C. - Rev.: 03/2012



RAPPRESENTANTE DI ZONA:
AREA SALES REPRESENTATIVE:



www.seedos.com

Medical Equipment and Services

