

**STERYLAB PRODUCT RANGE**

Biopsy & localisation needles  
Gynaecological products  
Decontaminating mats for the O.R.  
Aghi per biopsia e localizzazione  
Prodotti per ginecologia  
Tappeti decontaminanti

For further information or a quotation please contact : Colin Walters



**SeeDOS<sup>UK</sup>**

Radiation Oncology, Urology and General Medical Equipment and Services

26, The Maltings, Leighton Buzzard, Bedfordshire LU7 4BS, United Kingdom


Tel: +44 1525 850670 Fax: +44 1525 850685 Email: cwalters@seedos.com

Website: [www.seedos.com](http://www.seedos.com) Mobile: +44 07701059163

**PISTOLE PLURIUSO PER ASPIRATI CITOLOGICI**  
**REUSABLE SYRINGE HOLDERS FOR FINE NEEDLE ASPIRATION**

**PAC<sup>®</sup>**

■ Per uso con siringhe da 10, 20, 30 cm  
■ For use with 10, 20, 30 cm syringes




■ Per uso con siringhe da 10, 20 cm  
■ For use with 10, 20 cm syringes




**PAC-UF<sup>®</sup>**

# PISTOLE PLURIUSO PER ASPIRATI CITOLOGICI

## REUSABLE SYRINGE HOLDERS FOR FINE NEEDLE ASPIRATION

 Queste impugnature (pistole) consentono di effettuare aspirati citologici e svuotamenti di cisti con una sola mano, creando con facilità una notevole forza aspirante.

- ▶ Il modello PAC<sup>®</sup>, interamente costruito in resine plastiche e leghe leggere, consente l'utilizzo di siringhe da 10, 20 e 30ml delle marche più diffuse.
- ▶ Il modello PAC-UF<sup>®</sup>, costruito interamente in acciaio inossidabile e leghe leggere con specifici trattamenti superficiali che ne aumentano ulteriormente l'inattaccabilità da soluzioni disinfettanti aggressive, consente di utilizzare siringhe da 10 e 20 ml di qualsiasi marca.

 These useful accessories allow to perform cytologic sampling and emptying of cystis single handed, easily originating the vacuum into the syringe (not provided).

- ▶ Model PAC<sup>®</sup> is entirely made from plastics and light metal alloys and is designed to fit with the most common brands of 10, 20 and 30 ml syringes.
- ▶ Model PAC-UF<sup>®</sup>, which is entirely made from stainless steel and light metal alloys, with special surface treatment to increase resistance to even the most corrosive disinfectant, is designed to fit with the most common brands of 10 and 20 ml syringes.

### MODELLI DISPONIBILI

AVAILABLE TYPES

**PAC<sup>®</sup>**  
**PAC-UF<sup>®</sup>**

MODELLO MODEL	DESCRIZIONE DESCRIPTION	CODICE D' ORDINE ORDERING
PAC <sup>®</sup>	Modello per siringhe da 10, 20 e 30 ml model for 10, 20 and 30 ml syringes	PAC102030
PAC-UF <sup>®</sup>	Modello per siringhe da 10 e 20 ml model for 10 and 20 ml syringes	PACUF1020

Emiss.: D.C. - Rev.: 09/2004



RAPPRESENTANTE DI ZONA:  
AREA SALES REPRESENTATIVE:



[www.seedos.com](http://www.seedos.com)

Medical Equipment and Services

