


CENTERLINES REFER TO BOTH THE CENTROID OF THE COLLECTING VOLUME AND THE CENTER OF THE PHANTOM BLOCK

'F' BALL NOSE DRILL
x 5.985 DP



ALL DIMENSIONS IN INCHES AND TOLERANCES OF ± 0.005 UNLESS OTHERWISE NOTED

MATERIAL:
1. ACRYLIC, POLYSTYRENE, OR
OTHER; 30 cm x 30 cm

REVISION NO.

DRAWN BY:
BDH

PRODUCT
NOTE

SUBJECT : MODEL A14SL PHANTOM BLOCK MACHINING

DATE 11/1/01

PAGE 1 OF 1

NUMBER 4472-00

SeeDOS Ltd
cwalters@seedos.com
www.seedos.com